# CardioWatch Pro — ECG Analysis Pipeline Validation

On-device analysis of Apple Watch single-lead ECG recordings · zenmindstudio.com

## What the pipeline does

CardioWatch Pro processes 30-second Apple Watch ECG recordings (single lead, ~512 Hz) entirely on the user's iPhone: Butterworth band-pass filtering (0.5–40 Hz) → gradient-based R-peak detection → wavelet (DWT) wave delineation → interval measurement (RR, PR, QRS, QT; QTc via Bazett) → conservative rule-based labeling of signal features such as premature ventricular complexes (PVCs). Measurements are always displayed alongside clinical reference ranges. No data leaves the device.

## Validation datasets

- **MIT-BIH Arrhythmia Database** (PhysioNet) — records 100, 119, 207, 208, 233; thirty 30-second windows resampled from 360 Hz to 512 Hz to match Apple Watch sampling. Ground truth: cardiologist beat annotations.
- **MIT-BIH Normal Sinus Rhythm Database** — 5 records, thirty 30-second windows, for false-positive testing on clean rhythm.
- **76 real Apple Watch recordings** — spanning Apple's own classifications: sinus rhythm (60), high heart rate (8), poor recording (7), inconclusive (1).

## PVC labeling vs. cardiologist annotations

**79.4%**

Recall — share of annotated PVCs the pipeline caught

**94.8%**

Precision — share of PVC labels that matched an annotated PVC

**0.864**

F1 score (301 ground-truth PVC beats)

Beat matching tolerance:  $\pm 195$  ms against MIT-BIH beat-center annotations.

## False-positive behavior on normal rhythm

Recording set	Duration	False PVC labels	Rate
All normal-rhythm windows (MIT-BIH NSR + Apple Watch)	56 min	12	$\approx 1$ per 9 recordings
Apple Watch, classified "sinus rhythm" by Apple	30 min	9	$\approx 1$ per 7 recordings
Apple Watch, "high heart rate" / "poor recording"	7.5 min	0	0

One recording = one 30-second ECG window.

## Limitations

- Single-lead ECG cannot assess ST-segment changes; the pipeline deliberately performs no ST or ischemia analysis.
- Rule-based labeling is tuned to be conservative; recall is lower on records with unusual PVC morphology (MIT-BIH 207/208).
- Low signal-quality recordings are flagged as such rather than labeled; low-confidence results are marked "reference only" in the app.
- Validation is ongoing; this summary reflects the test suite shipped with app version 1.7.

**Important.** CardioWatch Pro is not a medical device and does not provide medical diagnosis. All measurements and labels are for reference only and cannot replace professional medical advice. If you experience chest pain, severe dizziness, or difficulty breathing, call emergency services immediately.

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